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*

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Read the package insert before use.
Retain the package insert and refer to it as needed.

In vitro diagnostic reagent
Identification number: 13A2X00044000002

** Revised in Aug. 2013(the ninth version)
* Revised in Feb. 2013(the eighth version)

Glycosylated hemoglobin A1c kit RAPIDIA® Auto HbA1c-L

General Precautions

1. This product is intended for in vitro diagnostic use only and not for other purposes.
2. Diagnosis should be made comprehensively on the basis of other relevant test results and clinical symptoms.
3. Any change or modification in the procedure described in the package insert may affect the reliability of test results.
4. Reagents in the Calibration Set for RAPIDIA® Auto HbA1c-L, sold separately, contain human-derived components. These are potentially infectious and should be handled accordingly.
5. Before use, read the package inserts and instruction manuals for devices that are used with these reagents.
6. This product contains sodium azide as a preservative. In case reagents accidentally contacted either with eyes, mouth or skin, clean up the contacted parts with plenty of water as a first-aid treatment, consult a doctor, if such is required.
7. HbA1c concentrations used in this package insert are expressed under the NGSP/JDS. The NGSP can be converted into the JDS using the following calculation formula.

$$\text{JDS}(\%) = 0.980 \times \text{NGSP}(\%) - 0.245\%$$

Kit Components

No.	Reagent	Form	Kit Size	
			20mL × 1	60mL × 1
R1	Latex	Liquid	60mL × 1	60mL × 3
R2	Anti-HbA1c Anti-human hemoglobin A1c mouse monoclonal antibody	Liquid	20mL × 1	60mL × 1
	Sample Diluent	Liquid	150mL × 1	150mL × 2

○ Product sold separately


Calibration Set for RAPIDIA® Auto HbA1c-L (5 concentrations × 3)

Intended Use

The test kit provides a quantitative assay for measuring concentrations of HbA1c in whole blood.

Principle of the Test

This kit is based on latex agglutination. HbA1c in test samples is absorbed onto the surface of latex particles, which react with Anti-HbA1c (antigen-antibody reaction). The turbidity caused by latex agglutination is measured at 660 nm, and the HbA1c concentration in whole blood is calculated from a standard curve.

Procedural Precautions

1. Properties and collection of test samples
 - 1) All samples should be collected by using EDTA-2K or NaF + Heparin + EDTA-2Na or NaF + EDTA-2Na contained collection tube.
 - 2) When measuring on the day of collecting blood samples, centrifuge the blood sample and collect packed red cell as a test sample. (700~1000G for 2 minutes)
 - 3) Blood samples should be stored at 2-10 °C and tested within two weeks. Do not freeze blood samples to avoid hemolysis.

2. Interfering substances and drugs

High lipid (3000 formazin turbidity units) or bilirubin (conjugated type 30 mg/dL, free type 30 mg/dL) counts in blood samples do not affect the test results.

3. Handling precautions

- 1) Use thoroughly washed test tubes to avoid errors in the assay.
- 2) Accurately dispense samples because the dispensing process influences the precision. Change tips each time to avoid errors due to cross-contamination.
- 3) Perform the assay step by step without delay.
- 4) The standard curve should be prepared for each set of assays since the reaction is affected to some extent by time and temperature, etc.

Usage and Dosage (assay procedures)

1. Preparation of reagents

No.	Reagent	Preparation	Stability
R1	Latex	Ready for use	1 month at 2 - 10°C after opening
R2	Anti-HbA1c	Ready for use	1 month at 2 - 10°C after opening
	Sample Diluent	Ready for use	1 month at 2 - 10°C after opening

*Mix each reagent well before use.

2. Materials required but not provided

- Pipette (1 mL)
- Micropipettes (2-20, 100, 300, 500 µL)
- Test tubes
- Mixer for test tube
- Autoanalyzer

3. Test procedures

For the assay, please refer to the operation manual of the autoanalyzer.

• Autoanalyzer: Example with a Hitachi 7170 (two-point end method)

1) Preparation of test samples (sample dilution)

The centrifuged blood sample, take 5 µL of packed red cells and hemolyze by adding 500 µL of Sample Diluent to prepare the test samples.

2) Solid phase reaction

Add 180 µL of R1 Latex to 4.8 µL of each test sample or each HbA1c calibrator. After mixing, incubate the mixture for 5 minutes at 37°C.

3) Antigen-antibody reaction

Add 60 µL of R2 Anti-HbA1c. After mixing, incubate for 5 minutes at 37°C.

4) Measurement

Measure the absorbance of test samples and respective HbA1c calibrators (main wavelength, 660 nm; secondary wavelength, 800 nm) according to the operating procedures for the autoanalyzer.

5) Calculation of concentrations

According to the operating procedures for the autoanalyzer, calculate the concentration of HbA1c in the test samples.

*When concentrations exceed the measurable range, dilute the test samples two-to three-fold with diluent (HbA1c calibrator L or another test sample of known concentration (%) of HbA1c). Run the test again and calculate the concentrations using the following formula:


$$\text{HbA1c concentration}(\%) = \frac{[\text{Retested HbA1c}(\%) \times n] - [\text{HbA1c}(\%) \text{ in the diluent} \times (n - 1)]}{n - 1}$$
 n = Dilution factor of the test sample (2 to 3)

[Example for calculation of concentrations exceeding the measurable range]

After the initial measurement, 100 µL of the test sample (HbA1c, approximately 17%) was diluted three-fold with 200 µL of a sample of known HbA1c concentration (4.8%) and retested, resulting in a value of 9.1%. In this case, the HbA1c concentration (%) is calculated as follows:

$$\begin{aligned} \text{HbA1c}(\%) &= \frac{(9.1 \times 3) - [4.8 \times (3 - 1)]}{3 - 1} \\ &= 27.3 - 9.6 \\ &= 17.7\% \end{aligned}$$

• Retested HbA1c (%) = 9.1%
 • n = 3
 • HbA1c(%) in the diluent = 4.8%
 • n - 1 = 2

Reference Normal values

HbA1c (NGSP) 4.6 – 6.2%

(HbA1c (JDS) 4.3 – 5.8%)

Performances **

1. Performance

1. Sensitivity

(1) When absorbance of Sample Diluent is measured according to the test procedure described in the package insert, the absorbance (turbidity) at 660 nm should be 0.005 or less.

(2) When absorbance of HbA1c calibrator H^{※1)} is measured according to the test procedure described in the package insert, the absorbance (turbidity) at 660 nm should be 0.15 or more.

2. Specificity

When three HbA1c control samples with a known concentration (%) (high, middle, low) are tested according to the test procedure described in the package insert, each test result should be within $\pm 10\%$ of respective known concentrations (%) of the HbA1c control samples.

3. Reproducibility

When 3 positive samples with a known concentration (high, middle, low) are tested times respectively according to the test procedure described in the package insert, the respective CV should be less than 10%.

4. Measuring range

HbA1c (NGSP) 3.3 ~ 13.5 %

(HbA1c (JDS) 3.0 – 13.0%)^{※2)}

※1 : HbA1c concentration of HbA1c calibrator H in Calibrator set for RAPIDIA[®] Auto HbA1c-L is 12% or more.

※2 : Result by Hitachi 7170.


2. Correlation

Study results in comparison with another latex agglutination method (in-house data)
n = 110

Correlation coefficient: $r = 0.999$

Regression equation: $y = 0.980x + 0.148$

(x, Approved latex method; y, RAPIDIA[®] Auto HbA1c-L)

3. Reference material for calibration

JCCRM411

Precautions for Use or Handling

1. Handling precautions (safety warnings)

- 1) Test samples should be handled with care, considering possible infection by HBV, HIV or HCV and others.

- 2) As in the case of samples, carefully handle calibrators in the Calibration Set for RAPIDIA[®] Auto HbA1c-L (sold separately), although the calibrators have been tested and found negative for HBs antigen, HIV antibody, and HCV antibody.
- 3) Latex and Sample Diluent in this kit contain sodium azide as a preservative. In case reagents accidentally contacted either with eyes, mouth or skin, clean up the contacted parts with plenty of water as a first-aid treatment, consult a doctor, if such is required.
- 4) Avoid freezing or drying the Latex because it may cause nonspecific agglutination.

2. Precautions for use

- 1) Do not use the reagents in this kit for purposes other than measurement of HbA1c in whole blood.
- 2) Strictly follow the storage conditions for this kit.
- 3) Do not use the reagents after the expiration date.
- 4) The reagents in this kit are arranged to ensure accurate reaction. Do not interchange reagents between different lot numbers.

3. Precautions for disposal

- 1) Samples, liquid waste, and equipment should be disinfected with hypochlorite or autoclaved at 121°C for at least 20 minutes.
- 2) Latex and Sample Diluent in this kit contain sodium azide (0.01% w/w) as a preservative. When disposing of the reagents, flush with large amounts of water to prevent the formation of explosive metal azide.
- 3) When disposing of the containers, reagents, liquid waste, and equipment, follow the applicable laws and regulations in each country.

Storage and Shelf Life

1. Storage

Store at 2 – 10°C .

2. Shelf life

Twelve months (refer to the expiration date printed on the carton and containers)

Package

1 kit: 20 mL × 1

60 mL × 1

○ Product sold separately

Calibration Set for RAPIDIA[®] Auto HbA1c-L (5 concentrations × 3)

Additional Information *

For additional information and further details consult:

Distributor: FUJIREBIO INC.

2-1-1 Nishishinjuku, Shinjuku-ku, Tokyo 163-0410 Japan

TEL: +81-3-6279-0899

Distributed by FUJIREBIO INC.

Manufactured by TFB, INC.